
Determination of Maximum Wavelength of Cr^{+6} Complex With Diphenylcarbazide (DPC)



Basic Principle:

It is little Violet (UV) worldle

absorption spectrophotometry.

Chemicals Required:

OIM Cott solution, 0.01% DPC

solution and 0.2N H2504


Observations and Colculations:

130	Wavelength	Absorbance	#ão	vovelength	Absentate
5.	450 nm	0.4	2	530 nm	1.71
2	460 nm	0.56	10	540 nm	1.83
3.	470 rm	0.78	11.	550 nm	1.80
4	480 nm	0 89	12	560 nm	1.75
5	490 11	1.09	13.	570 4	1.68
6.	500 nm	1 11	14.	580 mm	1 57
7.	510 nm	1.43	15.	590 **	1.49
8	520 mm	1.62	16.	600 nm	134

Result:

The Arms Value of Cx" DPC complex to

Spectroscopic Determination of Cr⁺⁶ Ions By Using Diphenylcarbazide in Water Sample



Basic Principle: It is UV visible absorption spectrscopu Chemicals: 0.1 M stock solution of Cot ions, 11. DPC solution. Acetone. Acetic acid/and distilled water. Preparation of Stock Solution: M = mass in grams x 1
molor weight volume of solution in dm3 $0.1 = \frac{\chi}{294} \times \frac{1}{0.275} \times = 0.275 \times 2.94 \times 0.1 = 8.085 \text{ gm}$ 8.085 gm K2 C62 O7 + Water up to 275 ml = 0.1 M K16,07 Preparation of Standard: MIVI = MIVI 0.1 x V1 = 0.005 X 100 $V_3 = \frac{0.005 \times 100}{0.1} = 5 \text{ m}$ 5 ml stock solution + water upto 100 ml = 0.005 M 10 ml " " + " = 0.01 M15 ml + + o = 0.015 M20 ml " " = 0.02 M= 0.025 M30 ml + + = 0.03 M

35 m	1 stock solution + Water	upto 100 ml = 0.035				
STATE OF THE STATE OF	, t					
1 2 2 4		- = 0.045 M				
50 ml	50 ml . + = 0.05 M					
Preparation of DPC Solution:						
1 gm DPC + 100 m/ acetone + CH3COOH						
Observations:						
S. No	Cx* Concentration	Absorbance				
1.	0.005 M	0.78				
2.	0.01 M	0.88				
3.	0.015 M	0.98				
4.	0.02 M	1.08				
5.	0.025 M	1.18				
6.	0 0 3 M	1.28				
7	0.035 M	1.39				
8.	0.04 M	1.47				
9.	0.045/M	1.58				
10.	0.05 M	1.7				
11	Sample (sud)	1.06				
Result: 28-11-16						
The concentration of Cx 10 ions in the						
given solution (sample) is 0.02 M.						